

SelectAir®



A complete line of
advanced & affordable
LOW AIR LOSS
THERAPY SYSTEMS

Featuring Two Computerized Control Units:

SelectAir® and SelectAir® "Max"

For use with our:

Low Air Loss Mattress Replacements
Low Air Loss Mattress Overlays

AirCare
Therapy Inc.

"Specialty Support Surfaces"



SelectAir® and SelectAir® “Max”

AirCare Therapy proudly introduces a complete line of Low Air Loss Therapy Systems. These unique, pressure-relieving support systems are designed to provide technologically advanced, computer-controlled, user-friendly operation. Portable and easily set up, our systems offer the most cost-effective alternative to pressure-sore management and prevention currently available.



SelectAir® Control Unit



SelectAir® “Max” Control Unit



SelectAir® and SelectAir® “Max” Mattress System

SelectAir® Cover Sheets are waterproof, breathable and designed to eliminate friction and shearing.

Low Air Loss Foot Pillow adds to patient comfort and prevents unnecessary skin breakdown due to patient migration.

33 Air Openings in the mattress promote the wicking of moisture.

SelectAir® Mattresses attach to any standard medical bed frame.

SelectAir®

- Reliable, Whisper-Quiet Operation
- Portable
- Lightweight (9 lbs.)
- Low Energy Consumption
- User-Friendly Control Panel
- Internal Usage Timer
- Microprocessor Controlled

Firm Mattress Mode instantly firms the mattress surface to facilitate nursing care during dressing changes and patient repositioning.

Upright Patient Mode gives added support when higher elevation is required, with no loss of clinical effectiveness.

SelectAir® “Max”

All the benefits of the *SelectAir* with these added features:

System Lock prevents unwanted adjustments, eliminating unnecessary service calls. Automatically resets the system in case of power failure.

During set up, **Auto Comfort** automatically sets optimal pressure below capillary closure based on the patient's weight.

Pulse Mode provides three levels of therapeutic, alternating pressure relief.

- The *SelectAir* and *SelectAir* “Max” are compatible with both our standard mattress replacement and our mattress overlay.
- Both the *SelectAir* and *SelectAir* “Max” provide superior pressure-sore management in a variety of health care settings.

Formerly AirCare Therapy Inc., the *SelectAir* line of products is now manufactured and distributed exclusively by The ROHO Group, Inc.

To place orders,
call The ROHO Group direct
1-800-851-3449
www.therohogroup.com

