

Jay Dermafloat LAL 1100

Low Air Loss Therapy



The Jay Dermafloat LAL 1100 is designed to work hard for you and your patients in all types of settings. A deluxe comforter and high-output control unit provide therapeutic benefit to patients suffering from or at risk of developing pressure ulcers. The LAL 1100 is one of the lightest, most durable and easiest to use genuine low air loss systems on the market.

Control Unit

- Easy-to-use dial adjusts pressure level from soft to firm.
- Compact, lightweight and durable case.
- Quiet operation allows patient to sleep more soundly.
- Blower-driven system provides “genuine low air loss” with all the benefits of significant air movement.
- CPR quick-disconnect for rapid deflation.

Standard Low Air Loss Mattress

- Max inflates within seconds.
- To reduce skin maceration, “genuine low air loss” is provided throughout all twenty air cells of support.
- Air cushions are constructed of high quality, low shear nylon to reduce skin breakdown.
- Convoluted foam base provides additional safety and comfort.
- Low-friction, vapor-permeable comforter provides a barrier to liquids, air and bacteria.
- Rolls into its own carrying case for convenience and portability.

Homecare Low Air Loss Mattress

- Orthoderm Elite Homecare Mattress is compatible with the new LAL 1000 Series control unit. Wide mattress also available.



Specifications

Description	Model	Dimensions:	Control Unit: 9" x 10" x 5" (8.5 lbs.) Standard Mattress: 80" x 34" x 10" Homecare Mattress: 80" x 35" x 10" Wide Mattress 83" x 39" x 10"
10" Standard System	360200	Electrical Rating: Agency Approvals: Warranty:	115/230V
Replacement Standard Mattress & Comforter	360240		UL 544/CE
10" Homecare System	360260	Maximum User Weight Capacity:	Controller - 2 year limited
Replacement Homecare Mattress & Comforter	350511		Mattress - 1 year limited
10" Wide Mattress System	360280		
Replacement Wide Mattress & Comforter	350611		
Control Unit	360220		300 lbs.